


## Blood Pressure Meter Design Using Microchip's PIC24F Microcontroller and Analog Devices

Author: Zhang Feng,  
Microchip Technology Inc.

### INTRODUCTION

A Blood Pressure Meter (BPM) is a non-invasive device used to measure blood pressure. This application note demonstrates the implementation of a digital blood pressure meter using Microchip's PIC24FJ128GC010 microcontroller and MCP6N11 instrumentation amplifier.

**FIGURE 1: BLOCK DIAGRAM**



## ANALOG FRONT-END SIGNAL CONDITIONING CIRCUIT DESCRIPTION

As shown in [Figure 1](#), the MCP6N11-100 single instrumentation amplifier (INA) is used to condition the analog signal acquired from the Wheatstone bridge type of air pressure sensor. The overall INA gain is set at 151 V/V. The output signal of the MCP6N11 is split into two paths. One path, representing the cuff pressure, is connected to the ADC1 channel of the microcontroller. Another path is passed through a 2-pole active high-pass filter with cutoff frequencies of 0.48Hz and 4.8Hz, and a gain of 92. The high-pass filter is formed by one of the microcontroller's internal operational amplifiers. The output of the high-pass filter, representing the oscillation signal, is sent to the ADC2 channel of the microcontroller through a notch filter.

## PIC24FJ128GC010 MICROCONTROLLER

The microcontroller provides a 12-bit high-speed pipeline Analog-to-Digital Converter (ADC). The cuff pressure signal and the oscillation signal are both sampled at 500 Hz by the 12-bit ADC.

As described in the previous section, the demo system uses one of the microcontroller's internal operational amplifiers to build the analog high-pass filter. It also utilizes one of the internal 10-bit Digital-to-Analog Converters (DAC) with buffered output voltage to provide an adjustable DC offset level as the bias for the operational amplifier.

Microchip's mTouch™ sensing solution is implemented in the firmware for the capacitive touch pads utilizing the microcontroller's Charge Time Measurement Unit (CTMU). Meanwhile, the air pressure sensor requires a stable supply current of 100 uA DC which is also provided by the CTMU.

The demo system utilizes the microcontroller's Pulse-Width Modulation (PWM) module to drive the air pump's motor in fast or slow speed mode. The noise generated by a running motor can interfere with the pressure measurement while inflating. The circuit design, firmware and Printed Circuit Board (PCB) layout are intended to reduce the motor noise.

The results for blood pressure and pulse rate can be sent to a Liquid Crystal Display (LCD) panel, USB interface and wireless interface. The microcontroller has an integrated LCD controller that generates the data and timing control required to directly drive a static or multiplexed LCD panel. The demo system utilizes the microcontroller's USB On-The-Go interface for data communication as well as the 5V main power source. For wireless communication, the microcontroller's SPI or UART port connects to one of Microchip's

Bluetooth® module RN42 or Wi-Fi® module RN171. A real-time clock is also displayed on the LCD panel by using the microcontroller's Real-Time Clock and Calendar module that can operate in Deep Sleep mode.

## POWER SUPPLY

The main input power is supplied by either four 1.5V AAA alkaline batteries or the USB's VBus line. The MCP1802 LDO converts the 5V or 6V main input voltage to 3.3V VDD.

## SPECIAL FEATURES

A digital low-pass filter with a cutoff frequency of 5 Hz is implemented in the firmware. It filters the raw ADC samples acquired on the oscillation signal.


This demo system uses a custom algorithm for blood pressure calculation that was developed based upon published academic papers available in the public domain (see [Appendix D: "References"](#)).

## MEASUREMENT TECHNOLOGY


A digital blood pressure meter measures systolic and diastolic pressures by oscillometric detection. Microchip's digital blood pressure meter demo can measure blood pressure and pulse rate during inflation. The Measurement While Inflating (MWI) principle reduces overall measuring time, which in turn reduces discomfort caused by the pressure in the cuff.

After the motor pumps the pressure up to 30 mmHg in fast mode, the motor changes to slow mode gradually and linearly inflating the cuff. When the pressure in the cuff gets high enough to just begin constricting the flow of blood, the arterial pulse becomes detectable by the pressure sensor. At this point, the waveform captured from the ADC2 channel (the output of the analog high-pass filter) starts to show the onset of the blood pressure oscillation signal (see [Figure 2](#)). The oscillation signal is then filtered by a digital low-pass filter. The Mean Arterial Pressure (MAP) is the cuff pressure which corresponds to the maximum oscillation signal. Based on the custom algorithm for blood pressure calculation, the Systolic Pressure (SYS) and the Diastolic Pressure (DIA) can be determined using the MAP value and a lookup table of empirical formulas (see [Figure 3](#)). The inflation process is automatically terminated when the pressure reaches a specified value. The Pulse Rate (PR) can be calculated based on ADC sample rate and sample numbers obtained in multiple consecutive pulses.

**FIGURE 2: RAW ADC DATA OF THE OSCILLATION SIGNAL ACQUIRED FROM THE ADC2 AT 500 SAMPLES/SECOND**




**FIGURE 3: THE PLOT OF ADC DATA FOR DIGITALLY FILTERED OSCILLATION SIGNAL AND CUFF PRESSURE SIGNAL**



## PROGRAM FLOWCHART

Figure 4 shows the firmware process flow.


**FIGURE 4: FIRMWARE PROCESS FLOWCHART**



## APPENDIX A: BPM DEMO BOARD

This appendix shows the BPM demo board.


**FIGURE 5: BPM DEMO BOARD**



## APPENDIX B: SCHEMATICS


This appendix shows the Blood Pressure Meter schematics.

### SHEET 1: MICROCHIP BLOOD PRESSURE METER SCHEMATIC – ANALOG FRONT-END




# SHEET 2: MICROCHIP BLOOD PRESSURE METER SCHEMATIC – PIC24F

## MICROCONTROLLER






## SHEET 3: MICROCHIP BLOOD PRESSURE METER SCHEMATIC – INTERFACE CIRCUITS







## **APPENDIX C: WARNINGS, RESTRICTIONS AND DISCLAIMER**

This demo is intended solely for evaluation and development purposes. It is not intended for medical diagnostic use.

## **APPENDIX D: REFERENCES**

Wang Wei-Wei, Pu Bao-Ming, He Bao-Yue, Li Sheng-Jin. Design and Realization of Electronic Sphygmomanometer Based on Inflation Measurement. Journal of Computer Systems & Applications, 2012,21(2):150-153,97.

---

**Note the following details of the code protection feature on Microchip devices:**

- Microchip products meet the specification contained in their particular Microchip Data Sheet.
- Microchip believes that its family of products is one of the most secure families of its kind on the market today, when used in the intended manner and under normal conditions.
- There are dishonest and possibly illegal methods used to breach the code protection feature. All of these methods, to our knowledge, require using the Microchip products in a manner outside the operating specifications contained in Microchip's Data Sheets. Most likely, the person doing so is engaged in theft of intellectual property.
- Microchip is willing to work with the customer who is concerned about the integrity of their code.
- Neither Microchip nor any other semiconductor manufacturer can guarantee the security of their code. Code protection does not mean that we are guaranteeing the product as "unbreakable."

Code protection is constantly evolving. We at Microchip are committed to continuously improving the code protection features of our products. Attempts to break Microchip's code protection feature may be a violation of the Digital Millennium Copyright Act. If such acts allow unauthorized access to your software or other copyrighted work, you may have a right to sue for relief under that Act.

---

Information contained in this publication regarding device applications and the like is provided only for your convenience and may be superseded by updates. It is your responsibility to ensure that your application meets with your specifications. MICROCHIP MAKES NO REPRESENTATIONS OR WARRANTIES OF ANY KIND WHETHER EXPRESS OR IMPLIED, WRITTEN OR ORAL, STATUTORY OR OTHERWISE, RELATED TO THE INFORMATION, INCLUDING BUT NOT LIMITED TO ITS CONDITION, QUALITY, PERFORMANCE, MERCHANTABILITY OR FITNESS FOR PURPOSE. Microchip disclaims all liability arising from this information and its use. Use of Microchip devices in life support and/or safety applications is entirely at the buyer's risk, and the buyer agrees to defend, indemnify and hold harmless Microchip from any and all damages, claims, suits, or expenses resulting from such use. No licenses are conveyed, implicitly or otherwise, under any Microchip intellectual property rights.

#### **Trademarks**

The Microchip name and logo, the Microchip logo, dsPIC, FlashFlex, KEELOQ, KEELOQ logo, MPLAB, PIC, PICmicro, PICSTART, PIC<sup>32</sup> logo, rPIC, SST, SST Logo, SuperFlash and UNI/O are registered trademarks of Microchip Technology Incorporated in the U.S.A. and other countries.

FilterLab, Hampshire, HI-TECH C, Linear Active Thermistor, MTP, SEEVAL and The Embedded Control Solutions Company are registered trademarks of Microchip Technology Incorporated in the U.S.A.

Silicon Storage Technology is a registered trademark of Microchip Technology Inc. in other countries.


Analog-for-the-Digital Age, Application Maestro, BodyCom, chipKIT, chipKIT logo, CodeGuard, dsPICDEM, dsPICDEM.net, dsPICworks, dsSPEAK, ECAN, ECONOMONITOR, FanSense, HI-TIDE, In-Circuit Serial Programming, ICSP, Mindi, MiWi, MPASM, MPF, MPLAB Certified logo, MPLIB, MPLINK, mTouch, Omniscient Code Generation, PICC, PICC-18, PICDEM, PICDEM.net, PICKit, PICTail, REAL ICE, rLAB, Select Mode, SQL, Serial Quad I/O, Total Endurance, TSHARC, UniWinDriver, WiperLock, ZENA and Z-Scale are trademarks of Microchip Technology Incorporated in the U.S.A. and other countries.

SQTP is a service mark of Microchip Technology Incorporated in the U.S.A.

GestIC and ULPP are registered trademarks of Microchip Technology Germany II GmbH & Co. KG, a subsidiary of Microchip Technology Inc., in other countries.

All other trademarks mentioned herein are property of their respective companies.

© 2013, Microchip Technology Incorporated, Printed in the U.S.A., All Rights Reserved.

 Printed on recycled paper.

ISBN: 978-1-62077-405-2

**QUALITY MANAGEMENT SYSTEM**  
**CERTIFIED BY DNV**  
**= ISO/TS 16949 =**

*Microchip received ISO/TS-16949:2009 certification for its worldwide headquarters, design and wafer fabrication facilities in Chandler and Tempe, Arizona; Gresham, Oregon and design centers in California and India. The Company's quality system processes and procedures are for its PIC® MCUs and dsPIC® DSCs, KEELOQ® code hopping devices, Serial EEPROMs, microperipherals, nonvolatile memory and analog products. In addition, Microchip's quality system for the design and manufacture of development systems is ISO 9001:2000 certified.*

## Worldwide Sales and Service

### AMERICAS

**Corporate Office**  
2355 West Chandler Blvd.  
Chandler, AZ 85224-6199  
Tel: 480-792-7200  
Fax: 480-792-7277  
Technical Support:  
<http://www.microchip.com/support>  
Web Address:  
[www.microchip.com](http://www.microchip.com)

**Atlanta**  
Duluth, GA  
Tel: 678-957-9614  
Fax: 678-957-1455

**Boston**  
Westborough, MA  
Tel: 774-760-0087  
Fax: 774-760-0088

**Chicago**  
Itasca, IL  
Tel: 630-285-0071  
Fax: 630-285-0075

**Cleveland**  
Independence, OH  
Tel: 216-447-0464  
Fax: 216-447-0643

**Dallas**  
Addison, TX  
Tel: 972-818-7423  
Fax: 972-818-2924

**Detroit**  
Farmington Hills, MI  
Tel: 248-538-2250  
Fax: 248-538-2260

**Indianapolis**  
Noblesville, IN  
Tel: 317-773-8323  
Fax: 317-773-5453

**Los Angeles**  
Mission Viejo, CA  
Tel: 949-462-9523  
Fax: 949-462-9608

**Santa Clara**  
Santa Clara, CA  
Tel: 408-961-6444  
Fax: 408-961-6445

**Toronto**  
Mississauga, Ontario,  
Canada  
Tel: 905-673-0699  
Fax: 905-673-6509

### ASIA/PACIFIC

**Asia Pacific Office**  
Suites 3707-14, 37th Floor  
Tower 6, The Gateway  
Harbour City, Kowloon  
Hong Kong  
Tel: 852-2401-1200  
Fax: 852-2401-3431

**Australia - Sydney**  
Tel: 61-2-9868-6733  
Fax: 61-2-9868-6755

**China - Beijing**  
Tel: 86-10-8569-7000  
Fax: 86-10-8528-2104

**China - Chengdu**  
Tel: 86-28-8665-5511  
Fax: 86-28-8665-7889

**China - Chongqing**  
Tel: 86-23-8980-9588  
Fax: 86-23-8980-9500

**China - Hangzhou**  
Tel: 86-571-2819-3187  
Fax: 86-571-2819-3189

**China - Hong Kong SAR**  
Tel: 852-2943-5100  
Fax: 852-2401-3431

**China - Nanjing**  
Tel: 86-25-8473-2460  
Fax: 86-25-8473-2470

**China - Qingdao**  
Tel: 86-532-8502-7355  
Fax: 86-532-8502-7205

**China - Shanghai**  
Tel: 86-21-5407-5533  
Fax: 86-21-5407-5066

**China - Shenyang**  
Tel: 86-24-2334-2829  
Fax: 86-24-2334-2393

**China - Shenzhen**  
Tel: 86-755-8864-2200  
Fax: 86-755-8203-1760

**China - Wuhan**  
Tel: 86-27-5980-5300  
Fax: 86-27-5980-5118

**China - Xian**  
Tel: 86-29-8833-7252  
Fax: 86-29-8833-7256

**China - Xiamen**  
Tel: 86-592-2388138  
Fax: 86-592-2388130

**China - Zhuhai**  
Tel: 86-756-3210040  
Fax: 86-756-3210049

### ASIA/PACIFIC

**India - Bangalore**  
Tel: 91-80-3090-4444  
Fax: 91-80-3090-4123

**India - New Delhi**  
Tel: 91-11-4160-8631  
Fax: 91-11-4160-8632

**India - Pune**  
Tel: 91-20-3019-1500

**Japan - Osaka**  
Tel: 81-6-6152-7160  
Fax: 81-6-6152-9310

**Japan - Tokyo**  
Tel: 81-3-6880-3770  
Fax: 81-3-6880-3771

**Korea - Daegu**  
Tel: 82-53-744-4301  
Fax: 82-53-744-4302

**Korea - Seoul**  
Tel: 82-2-554-7200  
Fax: 82-2-558-5932 or  
82-2-558-5934

**Malaysia - Kuala Lumpur**  
Tel: 60-3-6201-9857  
Fax: 60-3-6201-9859

**Malaysia - Penang**  
Tel: 60-4-227-8870  
Fax: 60-4-227-4068

**Philippines - Manila**  
Tel: 63-2-634-9065  
Fax: 63-2-634-9069

**Singapore**  
Tel: 65-6334-8870  
Fax: 65-6334-8850

**Taiwan - Hsin Chu**  
Tel: 886-3-5778-366  
Fax: 886-3-5770-955

**Taiwan - Kaohsiung**  
Tel: 886-7-213-7828  
Fax: 886-7-330-9305

**Taiwan - Taipei**  
Tel: 886-2-2508-8600  
Fax: 886-2-2508-0102

**Thailand - Bangkok**  
Tel: 66-2-694-1351  
Fax: 66-2-694-1350

### EUROPE

**Austria - Wels**  
Tel: 43-7242-2244-39  
Fax: 43-7242-2244-393

**Denmark - Copenhagen**  
Tel: 45-4450-2828  
Fax: 45-4485-2829

**France - Paris**  
Tel: 33-1-69-53-63-20  
Fax: 33-1-69-30-90-79

**Germany - Munich**  
Tel: 49-89-627-144-0  
Fax: 49-89-627-144-44

**Italy - Milan**  
Tel: 39-0331-742611  
Fax: 39-0331-466781

**Netherlands - Drunen**  
Tel: 31-416-690399  
Fax: 31-416-690340

**Spain - Madrid**  
Tel: 34-91-708-08-90  
Fax: 34-91-708-08-91

**UK - Wokingham**  
Tel: 44-118-921-5869  
Fax: 44-118-921-5820